


<b>Product Name :</b> Dynamics Of Machine Lab Equipments for Teaching Equipments Lab	<b>Product Code :</b> LABS-EDU-00138
	
<b>Description :</b>  <b>Dynamics Of Machine Lab Equipments</b>  <small>The set up consists of a small electric motor connected to a belt reduction gear box or belt and pulley system through a flexible coupling. The Governor spindle is driven by a belt coupling connected to the gear box output shaft and is supported on ball bearings. The governor mechanism can be mounted on a spindle. Speed control unit is provided to regulate speed and consequently the governor spindle can be used to run the mechanism. The governor spindle is in the gear and pulley mechanism. The weight down to which weights are added. The Governor can be operated as a static or portable Governor. The following range of experiments can be</small>	

Naugra Export

Website: www.naugraexport.com, Email: sales@naugraexport.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India, Phone: +91-0171-2643080, 2601773